

---

# Integrated Digital Microfluidic Functions for Chemical and Biological Applications

R.B. Fair

Department of Electrical and Computer  
Engineering

Duke University


Durham, N.C.



# Outline of Presentation


- Background and Motivation
  - Issues in true lab-on-a-chip systems
- Lab-on-a-Chip toolkit
  - Sample loading
  - Dispensing
  - Sample transport
  - Mixing reactors
  - Detection
- Performance data
  - Statistics on glucose assays
  - Pyrosequencing
- Summary and conclusions

# Background & Motivation



Test tubes

- ☐ Automation
- ☐ Integration
- ☐ Miniaturization




Robotics

- ☒ Automation
- ☒ Integration
- ☐ Miniaturization

Microfluidics

- ☒ Automation
- ☒ Integration
- ☒ Miniaturization




# Background and Motivation

- The reality of current lab-on-a-chip technologies...
  - Highly application specific
  - Commercial trend: simple, disposable devices that interface with expensive control boxes
  - Disposable devices may perform limited set of steps
- What is required for a true lab-on-a-chip?
  - Leverage devices into multiple applications
  - Complexity of diverse applications reduced to a manageable set of fluidic operations
  - Modular architecture gives flexibility of choosing fundamental operations
  - Top-down design




# Hierarchical Integrated Microfluidic Design


**Biomedical Fluidic Functions:**      **Func.1, Func.2,.....,Func.n**




**Elemental Set of Operations:**      **Op.1, Op.2,.....,Op.i**



**Elemental Set of Components**      **Comp. 1, Comp. 2,...,Comp. n**

- 
- Agent Detection
  - Precision Dispensing
  - Enzyme Analysis
  - Electrochromatography
  - Capillary Electrophoresis
  - Molecular/Protein Analysis
  - Isotachophoretic Separation

- 
- Transport
  - Mixing
  - Flushing
  - Filtering
  - Analysis
  - Detection
  - Monitoring

- 
- Buffers
  - Channels
  - Valves
  - Mixers

# Digital Microfluidic Toolkit

**Implementing numerous applications on a elemental set of components:**

Reservoirs → droplets

Dispensers → electrode sets

Pumps → electrode sets

Valves → electrode sets


Reaction vessels → droplets

Mixers → electrode sets

**Collection → scanning droplet**

# Implications of Droplet Architecture

- Droplets allow microfluidic functions to be reduced to a set of basic operations
- Numerous elemental fluidic operations can be accomplished with a common set of elemental components
- Array can be partitioned into “cells” that perform fluidic functions
- Functional cells dynamically reconfigured at least once per clock cycle




# Approach


- Develop a digital microfluidic lab-on-a-chip for analytical applications
- Specific focus
  - Multiplexed assays
  - Compatibility with biologically relevant liquids
  - Nanoliter scale
  - Sample-in-result-out operation
- Test lab-on-a-chip for clinical application
- Apply concepts to on-chip DNA sequencing



# Lab-on-a-chip Toolkit



**INTEGRATE**



Digital microfluidic  
lab-on-a-chip

# Sample Loading


**SAMPLE  
LOADING**

**DROPLET  
DISPENSING**


**DROPLET  
TRANSPORT**

**MIXING &  
REACTORS**

**DETECTION**



**TOP  
VIEW**



**SIDE  
VIEW**

- World-to-chip interface
- Loading using small volume pipette ( $<2\mu\text{L}$ )
- $W \ll R$  ensures that liquid stays in reservoir after loading

# Droplet Dispensing


**SAMPLE  
LOADING**

**DROPLET  
DISPENSING**

**DROPLET  
TRANSPORT**


**MIXING &  
REACTORS**

**DETECTION**



- High speed dispensing
- 8 droplets in 3.7 seconds
- Droplet volume ~ 20nL

# High Speed Continuous Droplet Dispensing





# Droplet Dispensing

**SAMPLE  
LOADING**

**DROPLET  
DISPENSING**

**DROPLET  
TRANSPORT**

**MIXING &  
REACTORS**

**DETECTION**

- Dispensing of proteins
  - Up to 1mg/mL BSA dispensable
  - Smaller than transportable concentrations
  - More adsorption due to larger surface area in the reservoir
- Dispensing of physiological fluids
  - Serum and plasma dispensable
  - Whole blood NOT dispensable
  - Use pressure-assisted dispensing

# Droplet Transport


**SAMPLE  
LOADING**

**DROPLET  
DISPENSING**

**DROPLET  
TRANSPORT**

**MIXING &  
REACTORS**

**DETECTION**



- High-speed transport
- 50Hz switching frequency
  - 2.5cm/sec speed
- 50V operation

# Droplet Mixing

SAMPLE  
LOADING

DROPLET  
DISPENSING

DROPLET  
TRANSPORT

MIXING &  
REACTORS

DETECTION

- Mixing in ~5 seconds by shuttling on linear array for 1  $\mu$ L (1.5mm scale) droplets
- Scaling down to 0.5mm will decrease mixing time
- Shuttling reverses flow causing un-mixing
  - unidirectional motion is preferred
- Mixing of two 25nL droplets was complete in 0.8 seconds at 10Hz switching @ 50V


# Rapid Droplet Mixing



2x4 Array Droplet Mixing

## Droplet Mixing on a 2x4 Electrode Array

Frequency:	16 Hz
Voltage:	50 V
Gap Height:	600 $\mu\text{m}$
Volume (each):	1.40 $\mu\text{l}$



- Droplets completely mix in 2.8 seconds
- 30 times faster than the diffusion-only passive mixing case



# Chemistry on Chip

**SAMPLE  
LOADING**

**DROPLET  
DISPENSING**


**DROPLET  
TRANSPORT**

**MIXING &  
REACTORS**

**DETECTION**

- Droplets are containerless chemical reactors
- Large and arbitrary dilutions are challenging and difficult to implement
  - Affects linearity of enzyme kinetic assays
  - Assays are sensitive to the sample composition
  - Possibility of interferences is much higher
  - Reactions requiring oxygen are affected
- Dilution factor of 2 (1 sample droplet + 1 reagent droplet) is most easily implemented


# Integrated Operation - Serial



- Serial protocol
- One glucose assay at a time
- Much simpler
- Does not require detector multiplexing

# Multiple Glucose Assays

Kinetic Data - Glucose Assays - 40, 80, 120 mg/dL




- 9 consecutive assays
- 3 glucose concentrations
- 60 seconds absorbance measurement at same spot

Data more noisy for 40 and 80mg/dL

# Multiple Glucose Assays

Serial Glucose Assays - within run variations



- *Potential* error sources: volume variation, cross contamination, measurement errors
- Low CV indicates good volume reproducibility
- No trend in the error indicates no cross contamination



# Detection Methodology




**SAMPLE  
LOADING**

**DROPLET  
DISPENSING**

**DROPLET  
TRANSPORT**

**MIXING &  
REACTORS**



**DETECTION**



# On-Chip Sequencing

- The promise of digital microfluidics to analyze large genome populations...
  - Miniaturized sequencing by synthesis using droplets
  - High throughput, massively parallel on chip reactions
  - Avoidance of by-product build-up that limits read lengths
  - Programmable operation
- Eliminating on-chip bottlenecks
  - Decouple the synthesis and detection steps
  - Introduce optimized concentrations of clean reagents for optimized times
  - Chemically amplified signal detection proceeds at leisure at a remote array site
- Potential to sequence an array of genomic segments each 10,000 to 100,000 bases long






STANFORD  
UNIVERSITY



Advanced Liquid Logic, Inc.



# Droplet-Based Pyrosequencing




Steps:

1. Move "A" nucleotide drop to DNA site and replace pyrobuffer drop
2. Move pyrobuffer drop to waste
3. Incubate incorporation reaction




# On-Chip Pyrosequencing




Steps:

1. Move incubated droplet and mix with luciferase
2. Move pyrobuffer drop and wash DNA (may be repeated)
3. Move combined droplet to array
4. Detect pyrophosphate generated light



# Example Fluidic Protocol for DNA Sequencing


- Dispense droplets of each dNTP
- Transport droplets to synthesis reaction site and allow to react
- Transport droplets to storage area
- Mix each dNTP droplet with light producing droplet
- Transport combined droplets to detector site





*Dyed liquids represent pyrosequencing reagents, droplet volumes are 50 – 100 nL.*

# On-chip Assay Development

On-chip pyrophosphate determination



dCTP incorporation into self-priming DNA substrate in solution



- Protocol: 1  $\mu$ L sample droplet mixed with 1  $\mu$ L pyrosequencing reagent droplet on-chip
- Pyrophosphate concentration is linearly correlated to luminescent peak height 0.3-3 pmol
- Self-primed DNA substrate + first complementary base (dCTP) is also linearly correlated
- Approximately 1 pmol DNA can be reliably detected with the current system

STANFORD  
UNIVERSITY



Advanced Liquid Logic, Inc.



# Lab-on-a-Chip Systems


- Digital microfluidic toolkit demonstrated
- Can digital microfluidics deliver a true lab-on-a-chip technology that is adaptable to numerous applications?
- Examples from ECE299 (Duke Univ. Fall 2006)
  - **Analog/Digital Hybrid Microfluidic Chip For DNA & RNA Analysis**
  - **Detection of DNA Using Fiber Optic Spectroscopy**






# Analog/digital hybrid biochip


(A. Garcia, G. Pan, J. Zhang)



# Fluidic Platform






# Floor Plan of the DMW



# Spectroscopic DNA Detection

(A. Vijay, N. Mathew, R. Erb)

- System that can detect the presence of specific DNA strands inside a blood sample in an entirely self-contained and hand-held system.




Aggregation is detectable!






# Fluid Level Flow Chart




## Block Diagram for Fluid Level



## Cell Lysing Area

In the Cell Lysis Area, Blood will be inputted in 25 nL sized droplets from the blood storage reservoir. This will be mixed with a 5 nL drop of rose buffer solution to lyse the cells. The end solution will be mixed with a 5 nL drop of magnetic beads. Then a 5 nL drop of solution will be pulled from the lysis area.



# Remarks on Applications

- Support concept that extensive biomedical application base can leveraged microfluidic operations in an electrowetting system.
- Based on:
  - Shared elemental fluidic operations
  - Reconfigurability
  - No cross-contamination
  - Multitasking by components
  - Few bottlenecks.
- Wide diversity of applications can be parsed into manageable components and assembled into a programmable, reconfigurable and reusable architecture.



# Summary and Conclusions

- Basic functionality of a true lab-on-a-chip demonstrated
  - Sample in/result out integration and automation achieved
- Electrowetting toolkit demonstrated
  - Automated droplet operations
  - Catalog of compatible reagents
  - Demonstrations of a few important biological assays
- Open issues:
  - On-chip sample preparation
  - System integration and interfacing to other laboratory formats and devices
    - Capillary electrophoresis
  - On-chip dilution still a problem
  - Scalable, compatible detector technology needed



# Acknowledgements

- NSF
- NIH
- DUHS
- ECE299 students

